Instructions for use

Disposable gloves
Supergloves Bluple nitrile, 868 Sizes S, M, L, and XL
Declaration of Conformity: www.carellurvink.nl/webshop/p/868h

Use

This instruction for use is to be used in combination with the specific information that is mentioned on or inside each packaging enclosure. These gloves are designed and should be disposed of after use. The gloves are liquid proof, and can therefore be used for splash protection against certain chemicals. They are in conformity with, and marked per the requirements of the European regulation 2016/425/EC and its amendments. They also comply with the applicable European Standards. Gloves which are accompanied with the pictogram which designates contact with foodstuffs are also in conformity with the European Regulations 1935/2004 as well as with applicable National Regulations for Food-contact materials. Please ensure the gloves are used only for the designated purposes.

Ingredients

Some gloves might contain ingredients which are known to be a possible allergies in sensitised persons, who may develop irritant and/or allergic contact reactions. If allergic reactions should occur, obtain medical advice immediately. For more information contact CLL B.V.

Care Instructions

Keep the gloves away from direct sunlight and store them in a cool, dry place. Keep the gloves away from ozone sources or a naked flame. The gloves are not designed to be laundered.

Disposal

Used gloves may be contaminated with infectious or other hazardous materials. Dispose of according to Local Authority Regulations. Landfill or incinerate under controlled conditions.

Precautions for use

- Before usage, inspect the gloves for any defects or imperfections. If the gloves are ripped or punctured during use, dispose of them immediately. If in doubt, do not use the gloves, get a new pair.

- If the gloves are used against chemicals, it is essential to keep all chemicals from the skin, even if they are thought to be harmless. Therefore use gloves which are rated with a protection index of 6 or with an excellent degradation resistance rating. In all other cases, the gloves should be used for splash protection or short contact only. Ensure that chemicals cannot enter via the cuff. For more details regarding chemical suitability please contact CLL B.V.

- Avoid wearing gloves which are dirty on the inside - they may irritate the skin causing dermatitis or worse.

- Contaminated gloves should be cleaned or washed before removal.

- The gloves should not be used in applications requiring thermal protection.

- Disposable glove type versions with a glove length below 260 mm are 'Fit For Special Purpose gloves' because they are to be protected the hand only from chemical splashes when handling chemicals. Do not use the gloves when protection in the cuffarea is needed.

- Gloves shall not be used for protection against ionising radiation nor for use in containment enclosures.

- Not all gloves that are suitable for contact with foodstuff can be used against all foodstuffs. Some gloves may show excessive migration towards certain types of foodstuffs. To know which restrictions apply and for which specific foodstuffs the gloves can be used, please obtain advise from CLL B.V.

- If gloves are marked, the printed surface shall not come in contact with food.

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This standard defines the general requirements for protective gloves in terms of construction, fitness of purpose, safety, etc.

Size s, m, l and xl
Some sizes comply with the EN 420: 2003 + A1: 2009 standard. Some gloves are shorter than a standard glove, to improve comfort for special purposes -

EU-Type examination certificate and Supervised product checks (Module C2) by
Centobel Gent
Technologiepark 7
9052 Gent
Tel. +32 9 220 41 51
Fax +32 9 220 49 55
gent@centobel.be

According to EN 420:2003+A1:2009, dexterity has been tested and reached level 5

Norm EN ISO 374-5 : 2016
Protective gloves against micro-organisms

The gloves must pass the permeation resistance test in accordance with the norm EN 374-2 : 2014. The possibility to offer resistance to viruses was added if a glove passed the following test, ISO 16604 : 2004 (method B).

EN ISO 374-5:2016

For gloves that offer protection against bacteria, fungi and viruses.

VIRUS

Norm EN ISO 374-1 : 2016
Protective gloves against dangerous chemicals

This norm is based on three testing methods

- * Permeation test in accordance with norm EN 374-2 : 2014
- * Permeation test in accordance with norm EN 16523-1 : 2015 which will replace norm EN 374-3
- * Degradability test in accordance with norm EN 374-4 : 2013

"This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals"

"The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture."

"It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation."

"When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves"

Six new chemicals were added to the list of harmful substances:

Code	Chemisch product	Code	Chemisch product
A	Methanol	K	40% natriumbisulfit
B	Aceton	L	90% zwavelzuur
C	Acetonitril	M	65% salpaterzuur
D	Dichloormethaan	N	99% salpzuur
E	Koolstofdioxide	O	25% ammoniak
F	Toluene	P	30% waterstofperoxide
G	Di-ethylamine	S	40% fluorzuur
H	Tetrahydrofuran	T	37% formaldehyde
I	Ethyl-acetaat		
J	n-Heptaan		

The following chemicals have been tested:

	Permeation Resistance EN 16523-1 ENISO 374-1:2016/type C	Determination of resistance to degradation by chemicals-perforation test(% degradation) EN ISO 374-4
K: NaOH, 40%	>=480 Level 6	3.2%
P: H2O2, 30%	>=15 Level 1	18.1%
T: formaldehyde, 37%	>=480 Level 6	39.6%

This standard defines the general requirements for protective gloves in terms of construction, fitness of purpose, safety, etc.

Size s, m, l and xl
Some sizes comply with the EN 420: 2003 + A1: 2009 standard. Some gloves are shorter than a standard glove, to improve comfort for special purposes -

EU-Type examination certificate and Supervised product checks (Module C2) by
Centobel Gent
Technologiepark 7
9052 Gent
Tel. +32 9 220 41 51
Fax +32 9 220 49 55
gent@centobel.be

According to EN 420:2003+A1:2009, dexterity has been tested and reached level 5

Norm EN ISO 374-5 : 2016
Protective gloves against micro-organisms

The gloves must pass the permeation resistance test in accordance with the norm EN 374-2 : 2014. The possibility to offer resistance to viruses was added if a glove passed the following test, ISO 16604 : 2004 (method B).

EN ISO 374-5:2016

For gloves that offer protection against bacteria, fungi and viruses.

VIRUS

Norm EN ISO 374-1 : 2016
Protective gloves against dangerous chemicals

This norm is based on three testing methods

- * Permeation test in accordance with norm EN 374-2 : 2014
- * Permeation test in accordance with norm EN 16523-1 : 2015 which will replace norm EN 374-3
- * Degradability test in accordance with norm EN 374-4 : 2013

"This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals"

"The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture."

"It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation."

"When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves"

Six new chemicals were added to the list of harmful substances:

Code	Chemisch product	Code	Chemisch product
A	Methanol	K	40% natriumbisulfit
B	Aceton	L	90% zwavelzuur
C	Acetonitril	M	65% salpaterzuur
D	Dichloormethaan	N	99% salpzuur
E	Koolstofdioxide	O	25% ammoniak
F	Toluene	P	30% waterstofperoxide
G	Di-ethylamine	S	40% fluorzuur
H	Tetrahydrofuran	T	37% formaldehyde
I	n-Heptaan		

The following chemicals have been tested:

	Permeation Resistance EN 16523-1 ENISO 374-1:2016/type C	Determination of resistance to degradation by chemicals-perforation test(% degradation) EN ISO 374-4
K: NaOH, 40%	>=480 Level 6	3.2%
P: H2O2, 30%	>=15 Level 1	18.1%
T: formaldehyde, 37%	>=480 Level 6	39.6%

CLL B.V.
IJzersteden 11
7547 TB Enschede

Tel: +31 (0) 53 434 43 43
Fax: +31 (0) 53 433 71 05
Mail: info@carellurvink.nl

KVK 67631266
BTW NL857100671B01
IBAN NL15ABNA0538346604
BIC ABNANL2A

CLL B.V.
IJzersteden 11
7547 TB Enschede

Tel: +31 (0) 53 434 43 43
Fax: +31 (0) 53 433 71 05
Mail: info@carellurvink.nl

KVK 67631266
BTW NL857100671B01
IBAN NL15ABNA0538346604
BIC ABNANL2A