





Easy to open and simple to use

- · Available with or without Filter, sized to fit any 400 ml Homogenizer machine
- Filter allows for separation of the liquid extract from sample solids, thereby preventing pipette blockage
- · Filter with solid debris maybe easier removed if desired
- · Bags slide easily for quick access

Strong enough to endure the sample Preparation process, even at low temperatures

- Thick 3 mil polyethylene for toughness
- Additional additive for increased elasticity
- Extra wide 4 mm heat seal is designed to prevent leaking
- Filter section of the bag is made of durable, perforated polythene to give extra reinforcement and support during the blending process



No worries about prior contamination of the bag

Gamma irradiated

Backed by 3M's quality reputation

• Manufactured with proven experience and uncompromising quality

For more product information, visit: www.3m.com/microbiology

Ordering Information

Product	Catalog No.	Contents	Dimensions
3M Filter Bag	6469	200 bags/case	180 x 300 mm
3M Homogenizer Bag Filtered 400	6468	500 bags/case	180 x 300 mm



3M Health Care Food Safety Department Europe, CEE, MEA Hammfelddamm 11 D - 41453 Neuss Tel. +49 (0)2131-14 3000 Fax: +49 (0)2131-14 4443 www.3m.com/microbiology